**IMPORTANT MEDICAL INFORMATION
SURGICAL INSTRUMENT**

Before using a product introduced on the market by Ortopro, the operating surgeon is asked to carefully study the following recommendations, warnings, and instructions

This user instruction applies for surgical instruments (hereinafter called "instrument") of Ortopro that are needed for the implantation of prostheses, for the use of joint replacement, internal & external fixation products.

1-) Safety Instructions

1.1 General Instructions

Products of Ortopro may only be used by surgeons who are conversant with the general problems of joint replacement, internal & external fixation surgery and who master the product-specific surgical techniques.

The instruments are to be stored either in their original packaging, unopened, in the appropriate instrument tray or in the suitable packing that protects the instruments from damage. Before they are used they must be examined for possible damage and it must be checked that they are in working order.

Any complications or other effects that may occur for reasons such as incorrect indication or surgical technique, unsuitable choice of material or treatment, inappropriate use or handling of the instruments etc., fall under the responsibility of the operating surgeon and can not be blamed on the manufacturer , the importers.

The Manufacturer does not accept liability in the case of non-compliance with the requirement of this user instruction.

Important instructions:

- Instruments of Ortopro are always components of a system. They may only be combined with original parts belonging to the same system.
- Instrument must not be mechanically worked or altered in any other way,
- The general risk involved in the use of instruments are allergic reaction to the material used in their manufacture, loosening, wear, corrosion, ageing and fracture of the instrument or its components.
- The use of the instruments for other purposes is forbidden.

Warnings:

- Instruments that are contaminated, unsterile or damaged or that have been improperly handled or instruments or their components that have been altered without authorization must not be used under any circumstances.
- If instruments are subjected to excessive loading, are damaged or are handled improperly, they may fracture, corrode, oxidize, become loose and show excessive wear or their functionality may be impaired.
- With the use of reamers, drills and other cutting instruments, frictional heat can be generated, which may lead to cell damage.
- Instruments may sometimes have sharp edges or pointed tips. Consequently, if plastic gloves are used there is a danger that these may become damaged. Therefore please pay attention to the risk of infection.
- Articulating or assembled instruments or their components must be thoroughly cleaned before they are assembled or used. Because contamination can lead to wear, impairment of function or fracture of the instruments or their instruments.
- Torque wrench accuracy is $\pm 10\%$ of set value over the life of the instrument, as long as it is recalibrated every 6 months by authorized Ortopro personnel. the torque wrench mechanism will relax, whether it is sitting on a shelf or being used in daily procedures. Ortopro therefore requires that each torque wrench be returned every six (6) months for a review of all components for wear and recalibration to maintain device accuracy.

1.2 Cleaning and Maintenance

Instruments must not come into contact with substances containing chlorine or fluorine. Instruments that are made wholly of plastic must not come into contact with strong acids or organic or ammonia-containing solvents, aromatic and/or halogen hydrocarbons or oxidising chemicals. Aluminium and materials containing aluminium must never come into contact with substances containing mercury. Even the smallest traces mercury can lead to considerable corrosion. Instruments made from materials containing aluminium must be wiped and cleaned with, or placed in detergents and disinfections that have a pH value between 4.5 and 8.5. with higher or lower pH values the protective neutral coating of materials containing aluminum dissolves, which leads to corrosion. Detergents containing caustic may cause the surface corrossions and discoloration on anodized aluminum Instrument set therefore Detergent containing caustic should not be used during cleaning process for those parts. Multiple-use instruments must be cleaned and sterilized after use (also see instructions under point 1.3)

Contaminations must not be allowed to dry on the instruments, otherwise, cleaning can be difficult.

Instrument should not be placed in physiological saline solution, as prolonged contact with this medium can lead to corrosion and changes to the surface of the instruments.

Freshly prepared cleaning and sterilizing materials must always be used. Neither metal brushes nor scouring agents may be used for the cleaning of instruments. In order to avoid water stains, a final rinsing with desalinated water is recommend. The instruments must then be dried immediately.

1.3 Sterilization

Sterilization must be carried out with equipment and under conditions that meet the requirements of applicable standards . the user is responsible for the regular maintenance and checking of the cleaning and sterilization equipment and for the appropriate validation.

The components are supplied non-sterile and must be sterilized prior to use by a hospital validated steam autoclaving process in appropriate protective wrapping when necessary. If necessary Components must be cleaned prior to sterilization in compliance with hospital validated cleaning process or Cleaning equipment manufacturers' user instructions and recommendations for chemical detergent is required. All instruments must be stored in a clean, dry environment and be protected from sunlight and extremes in temperature.

The following process parameters are validated by Ortopro and recommended for sterilization.


Steam Sterilization		
Cycle Type	Parameter	Minimum Set Point
Prevacuum 273°F (134°C)	Exposure Temperature	273°F (134°C)
	Exposure Time	4 minutes
	Dry Time	20 minutes

These recommendations have been developed and validated using specific equipment. Due to variations in environment and equipment, it must be demonstrated that these recommendations produce sterility in your environment.

Important instruction:

- **The validation of the cleaning, sterilization and reesterilization procedures and the correct setting of the corresponding equipment must be checked out regularly.**
- **The recommendations set out under points 1.2 and 1.3 serve only for information. No liability is accepted for instruments that are cleaned and sterilized by the purchase or the user, or for reesterilised instruments.**

2-) SYMBOLS

	Batch code		Date of manufacture
	Catalog number		contents packed without sterilization
	Consult Instruction For Use		

INFORMATION

Should any incident occur with implantable device, call the phone number given below. For further information, please contact Customer Service

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